
















Connected Care: Insulin Pump, CGM, and App Guide



AID System	CGM	Age and Diagnosis	Algorithm Target	Reservoir	Phone	Mobile App	Simulator
Omnipod Dash 		2+ T1D 18+ T2D	<u>No Algorithm</u> 110 - 150	 Max Insulin: 200 units Min. Insulin: 85 units Min. TDD: 5 units Min. Weight: N/A			
Omnipod 5 	  	2+ T1D 18+ T2D	<u>SmartAdjust</u> 110 - 150	 Max Insulin: 200 units Min. Insulin: 85 units Min. TDD: 5 units Min. Weight: N/A	 		
Tandem t:slim X2 	  	2+ T1D 18+ T2D	<u>Control IQ+</u> 112.5 - 160	 Max Insulin: 300 units Min. Insulin: 50 units Min. TDD: 5 units Min. Weight: 20 lbs./9 kg	 		
Tandem Mobi 	 	2+ T1D 18+ T2D	<u>Control IQ+</u> 112.5 - 160	 Max Insulin: 200 units Min. Insulin: 35 units Min. TDD: 5 units Min. Weight: 20 lbs./9 kg			

Connected Care: Insulin Pump, CGM, and App Guide



Insulin Pump	CGM	Age and Diagnosis	Algorithm Target	Reservoir	Phone	Mobile App	Simulator
iLet Bionic Pancreas 		6+ T1D	Go Bionic 110, 120, 130	 Max Insulin: 180 units (user-filled) 160 units (pre-filled*) Min. Insulin: N/A Min. TDD: N/A Min. Weight: N/A			
Sequel twist 		6+ T1D	Loop 87 - 180	 Max Insulin: 300 units Min. Insulin: 100 units Min. TDD: N/A Min. Weight: N/A			
MiniMed 780G 		7+ T1D 18+ T2D	SmartGuard 100 - 150	 Max Insulin: 300 units Min. Insulin: 100 units Min. TDD: N/A Min. Weight: N/A			

*iLet Bionic Pancreas pre-filled reservoir uses Fiasp (Ultra rapid-acting insulin aspart).



*Indicates Dexcom G7 10-day and 15-day sensor.

Pump/CGM Company	Contact Info
Insulet (Omnipod Dash; Omnipod 5)	800-591-3455 opt. 1
Tandem (t:slim X2; Mobi)	877-801-6901
Beta Bionics (iLet)	855-745-3800 opt. 1
Sequel MedTech (twiist)	877-489-4478
Medtronic (MiniMed 780G)	800-646-4633 opt. 1
Dexcom (G6; G7)	888-738-3646
Abbott (Freestyle Libre 2 Plus; 3 Plus)	855-632-8658